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Abbreviations

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References

1. Raman thermometry. **Tuschel, David**. Dec. 2016, Spectroscopy solutions for material analysis, pp. 8–13.
2. Raman temperature measurement. **D. S. Moore, S. D. McGrane**. 2014, J. Phys. Conf Ser., Vol. 500.
3. High temperatures and Raman scattering through pulsed spectroscopy and CCD detection. **P. Simon, B. Moulin, E. Buixaderas, N. Raimboux, E. Herault, B. Chazallon, H. Cattey, N. Magneron, J. Oswalt, D. Hocrelle**. 2003, J. Raman Spectrosc., Vol. 34, pp. 497–504.
4. Nanosecond time-resolved Raman spectroscopy for solving some Raman problems such as luminescence or thermal emission. **E. S. Fotso Guetue, A. Canizares, P. Simon, N. Raimboux, L. Hennet, M. -R. Ammar**. 2018, J. Raman Spectrosc., Vol. 49, pp. 822–829.
5. A 2x(4)x128 Multi-time-gated SPAD Line Detector for Pulsed Raman Spectroscopy. **I. Nissinen, J. Nissinen, P. Keränen, A. Lämsman, J. Holma, J. Kostamovaara**. 2015, IEEE J. Sens., Vol. 15, p. 1358.
6. On-line monitoring of spodumene heat treatment process with time-gated Raman spectroscopy. **P. Tanskanen, B. Heilala, L. Kurki, J. Savela**. 2018, Process Mineralogy '18 proceedings, Cape Town, South Africa.
7. Advantages of anti-Stokes Raman Scattering for high-temperature measurements. **Hirofumi Fujimori, Masato Kakihana, Koji Ioku, Seishi Gota, Masahiro Yoshimura**. 2001, Appl. Phys. Lett., Vol. 79, pp. 937–939.
8. Advantages of ultraviolet Raman scattering for high temperature investigations. **E. Zoubolis, D. Renusch, M. Grimsditch**. 1998, Appl. Phys. Lett., Vol. 72, pp. 1–3.
9. UV Raman spectroscopic characterization of catalysts and catalytic active sites. **Shaoqing Jin, Zhaochi Feng, Fengtao Fan, Can Li**. 2015, Catal. Lett., Vol. 145, pp. 468–481.
10. A comparison of ultraviolet and visible Raman spectra of supported metal oxide catalysts. **Yek Tann Chua, Peter C. Stair, Israel E. Wachs**. 2001, J. Phys. Chem. B, Vol. 105, pp. 8600–8601.
11. Temperature-dependent photodegradation in UV-resonance Raman spectroscopy. **Hikaru Yoshino, Yuika Saito, Yasuaki Kumamoto, Atushi Taguchi, Prabhat Verma, Satoshi Kawata**. 2015, Anal. Sci., Vol. 31, pp. 451–454.
12. High-temperature Raman spectroscopy. **Osipov, Armenak A.** 2019, Pure Appl. Chem., Vol. 91, pp. 1749–1756.
13. Raman spectroscopy at high pressures. **Goncharov, Alexander F.** 2012, Int. J. Spectrosc., Vol. 2012.

Author

This blog was written by Timegate Instruments' Senior Application Specialist **Bryan Heilala**. Bryan is a young and energetic chemist with a degree in M.Sc. (chemistry) and experience and background in analytical chemistry. Read more about him and the whole [Timegate team](#).